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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow its instructions may pose a threat to life or health.

APPLICATION AND DESCRIPTION OF THE DEVICE

This is a fully automatic digital blood pressure monitor designed for home use. It allows very fast and accurate measurement of systolic and diastolic blood pressure, as well as heart rate using the oscillometric measurement method. Our product does not use mercury or a stethoscope, which makes it safe and easy to use.

Please read this user manual carefully before use and keep it in a safe place. Consult your doctor if you have any questions about blood pressure and measurements.

For safety and CE certification reasons, the product may not be rebuilt or modified in any way. If the product is used for purposes other than those previously described, the product may be damaged. Improper use may also cause hazards such as short circuits, fire, electric shock, etc.

REMEMBER

- Only a qualified healthcare professional can interpret blood pressure measurement results.
- This device is NOT intended to replace regular medical examinations.
- It is recommended that your physician approve your use of the device.
- Blood pressure readings obtained with this device should be confirmed by a physician before changes are made to medications or dosage.
- Never change medication doses on your own based on measurements without consulting your doctor.
- This device is intended for adults only. Consult your doctor before using it on children.
- If arrhythmia (irregular heart rhythm) occurs, the results may be abnormal and should be evaluated by a physician.

Battery replacement instructions:

- Please insert the batteries correctly, paying attention to the +/- polarity.
- Only use one type of battery.
- Always replace the entire set of batteries.
- Exhausted batteries must be removed from the device and never thrown into fire.
- Never attempt to charge regular batteries.
- Store unused batteries in their packaging away from metal objects.
- Do not leave children unsupervised.

TECHNICAL DATA

- Display: LCD
- Pressure height: 70KPa~106KPa
- Measurement method: oscillometric
- Pressure sensor: capacitive
- Measurement range: SYS/DIA 0 to 280 mmHg / Pulse: 40 to 199 per minute
- Pressure display range: 0-299 mmHg
- Memory: up to 99 measurements
- Measurement resolution: 1mmHg
- Accuracy: Pressure: ± 3 mmHg / Pulse: $\pm 5\%$ max
- Automatic air pumping
- Automatic shutdown after 60 seconds
- Watch
- Date
- Simple operation
- Perfect for traveling
- Functional case
- Power supply: 4 AAA batteries (not included)
- Anti-slip base
- Dimensions of the device without cuff: 10.5cm x 12.5cm
- Display size: 6.5cm x 5.5cm
- Cuff dimensions: 48cm x 12.5cm

IMPORTANT INFORMATION ABOUT BLOOD PRESSURE AND ITS MEASUREMENT

Where does high or low blood pressure come from?

- Blood pressure levels are set by the central nervous system and adjust to different situations through feedback. To adjust blood pressure, the strength and speed of the heart (pulse) are changed, as well as the width of the blood vessels. The width of the blood vessels is controlled by tiny muscles in the walls of the blood vessels.
- Blood pressure levels change periodically during heart activity: During the "blood ejection" (systole), the value is the highest (systolic pressure value). At the end of the heart's "resting period" (diastole), the pressure is the lowest (diastolic pressure value).

Blood pressure values must be within certain ranges to prevent certain diseases.

X – Systolic blood pressure (mmHg)

Y – Diastolic blood pressure (mmHg)

A – Optimal blood pressure


B – Normal blood pressure

C – Increased blood pressure

D – Mild high blood pressure

E – Moderate blood pressure

F – Severe high blood pressure


If your blood pressure is normal, it is recommended that you regularly self-monitor your blood pressure using a blood pressure monitor.

You can detect any changes in blood pressure early and respond appropriately. If you are undergoing medical treatment to control your blood pressure, record the values along with the time of day and date.

Show these values to your doctor. Never change your doctor's prescribed medication without consulting your doctor.

DEVICE DESCRIPTION

BATTERY INSTALLATION

1. Remove the battery cover on the bottom of the device by sliding it in the direction shown on the cover.
2. Install 4 x “AAA” batteries so that the + (positive) and (negative) poles match the poles of the battery compartment. (Long-life alkaline batteries are recommended)

USING MAINS POWER

Warning: Do not install any battery when using the AC adapter.

1. Using the micro USB cable, connect the device to the AC adapter.
2. After taking measurements, disconnect the AC adapter.

FUNCTION USER MANUAL

1. Selecting a user

Press the SET button for at least 3 seconds, the symbol should flash on the display. Press the MEM button to select the user. Press the SET button to confirm.

2. Setting the time and date (year, month, day, hour and minute)

Perform part 1. Select a user. After selecting a user, the year on the display will flash, press the MEM button to increase the number by one (press and hold the MEM button to increase faster). to increase the number by one (press and hold the MEM button to increase faster). Press the SET button to set the current year, then the month on the display will flash.

Repeat the above procedure to set the month, day, hour and minute.

Note: There are 2 time modes 24-hour/12-hour, the initial mode is 24-hour, press the memory button to set the time mode, when you select 12-hour mode, the time will be displayed with AM and PM. Press the setting button to confirm.

3. Unit Setting: This monitor can switch the units between mmHg and kPa. After setting the time, PR OFF or PR ON will flash. Press the MEM button to switch. PR OFF indicates the unit is mmHg. PR ON indicates the unit is kPa. Press the SET button to confirm.

4. Connecting the cuff tube

Insert the cuff tube into the hole on the left side of the monitor, marked with a cuff drawing.

MEASUREMENT

1. Before measurement:

Do not eat, smoke or exercise before taking a measurement. These activities will affect the accuracy of the result. Relax for 10 minutes before taking a measurement and relax for 30 minutes after intense exercise.

Remove any clothing from your arm. Do not roll up your sleeve, this can increase the pressure on your arm and affect the result,

Always measure on the same arm (usually the left one, close to the heart)

2. Applying the cuff:

Place the cuff on your left arm so that the air tube is on the inside of your arm and is aligned with your middle finger. The lower end of the cuff should be approximately

2 cm above the elbow (about two finger widths)

Wrap the cuff tightly, allowing only one finger (approximately 1 cm wide) to fit between the cuff and the arm.

2-3cm 22-32cm

3. Measurement procedure:

Sit in a chair, back straight. Place your forearm flat on the table, palm facing up.

Press the START/STOP button to start the measurement. The cuff starts to inflate automatically.

The LCD display shows the increasing pressure level in the cuff. The monitor determines the ideal inflation level and stops when this level is reached. Then it starts to detect the pulse.

Once the measurements are taken, the cuff will begin to deflate.

The display shows blood pressure and heart rate.

4. Press the START/STOP button to turn off the monitor. It will turn off automatically after one minute.

Note: Do not move your body in any way until the entire measurement is completed.

MEMORY FUNCTION

At the end of the measurement, the device automatically saves the result with date and time. The monitor stores 99 readings for 2 users (user 1 and 2)

Results overview:

1. Select a user
2. Press the MEM button and the screen will display the last measurement result. The first reading displayed will be the average of the last 3 results. Press the MEM button again and the monitor will display the previous result from newest to oldest.
3. Once measurement is performed more than 99 times, reading No. 99 will be deleted and the new reading will be saved.

Memory deletion:

1. Delete a single reading:

When the reading you want to delete appears on the display, hold down the SET button until the upper right corner should read "DEL". Press the MEM button to delete the reading.

2. Delete all readings:

Hold the MEM button for at least 3 seconds (don't release), press the SET button, then the display will show a message. Stop holding the MEM and SET buttons at the same time, and all readings will be cleared.

INTERRUPTING MEASUREMENT

If for any reason it is necessary to interrupt the blood pressure measurement, the START/STOP button can be pressed at any time. The device will immediately and automatically lower the pressure in the cuff.

IRREGULAR HEARTBEAT INDICATOR

When the monitor detects an irregular rhythm, an irregular heartbeat symbol will appear in the lower left corner of the screen.

An irregular rhythm is a rhythm that differs by less than 25% from the average rhythm or by more than 25% from the average rhythm detected by the device when measuring systolic and diastolic blood pressure. If the irregular heartbeat symbol appears, consult your doctor.

ERRORS DESCRIPTION

Mistake	Cause	Solution
No display appears	Weak battery or improper battery placement	Replace both batteries with new ones. Check that the battery installation is correct.
Er P	The cuff does not inflate	Check for cuff leaks or insert air plug securely.
Er 1	Air escapes too quickly or the pulse is too weak	Check the cuff position and measure again.
Er 2	Signal interference	Remove anything that may interfere with heartbeat detection. (Cloth, electromagnetic wave, etc.)
Er 3	Incorrect result	Put the cuff on correctly and measure again
HI	Cuff over inflated	Correct the cuff position and repeat the measurement.

SAFETY INSTRUCTIONS

- This product may be used by children aged 8 years and over and by persons with reduced physical, sensory or mental capabilities or lack of experience or knowledge, provided they are supervised or have been instructed on the safe use of the product and understand the resulting risks.
- The device contains sensitive electronic components. Avoid strong electromagnetic fields (e.g. mobile phones, microwave ovens) that may cause interference.
- Do not use cuffs, adapters or batteries that do not meet technical requirements.
- Remove the battery if the device will not be used for a long time.
- Do not use near flammable gases.
- Cuff material is not suitable for children under 3 years of age.
- If different patients use the same device, wash your hands before and after each use.
- The cuff contains a sensitive, sealed bulb. Handle the cuff with care and avoid any type of stress by twisting or bending.
- Clean the device with a soft, dry cloth. Do not use gas, thinners or similar solvents. Stains on the cuff can be carefully removed with a damp cloth and soap suds.
- Handle the tube with care. Do not pull on it. Do not allow the tube to kink and keep it away from sharp edges.
- Avoid strong vibrations for the device.
- In case of arrhythmia – measurements should be evaluated by a physician.
- Dispose of the device according to local environmental regulations.

- Never immerse the device in water.
- Turn your device off after use to avoid unnecessary battery consumption.
- Check the batteries regularly and replace them when the device stops working properly. Dispose of used batteries according to local regulations.
- Use a soft, damp cloth for cleaning. Avoid using harsh cleaning agents that may damage the device.
- Do not use a damaged device.
- The product is not intended for use by children.
- Make sure that packaging materials are not left unattended. Children may start playing with them, which is dangerous.
- Protect the product from extreme temperatures, direct sunlight, strong vibrations, high humidity, moisture, flammable gases, vapors and solvents.
- Do not expose the product to mechanical stress.
- If safe operation is no longer possible, discontinue use and secure the product against further use. Safe operation is no longer possible if the product: - has been damaged, - does not function properly, - has been stored for an extended period under unfavourable conditions or - has been excessively stressed during transport.
- It is prohibited to use the product if any part is damaged.
- This product contains small parts that may present a choking hazard. Keep out of reach of children.
- Do not disassemble the device yourself.
- Do not use old and new batteries, different brands or chemicals.
- Do not charge, damage or open the battery.
- Keep batteries out of the reach of children to prevent accidental swallowing.
- Do not swallow the battery. If the battery is swallowed, consult a doctor immediately.
- Always use the product as intended.

BATTERIES/ACCUMULATORS SAFETY INSTRUCTIONS

- Batteries/rechargeable batteries should be kept out of reach of children. If swallowed, consult a doctor immediately!
- Disposable batteries must not be recharged. Batteries/accumulators must not be short-circuited and/or opened. This may cause overheating, fire or explosion.
- Never throw batteries/accumulators into fire or water.
- Batteries/accumulators should never be subjected to mechanical stress.
- Risk of acid leakage from batteries/accumulators.
- Avoid extreme conditions and temperatures that may affect batteries / accumulators, e.g. radiators / direct sunlight.
- If batteries/accumulators leak, avoid contact of skin, eyes and mucous membranes with chemicals! Immediately rinse affected areas with clean water and seek medical attention!
- Leaking or damaged batteries/rechargeable batteries can cause chemical burns when in contact with skin. Therefore, suitable protective gloves should be worn in such cases.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a waste disposal point designated by local authorities. Information on the possibilities of recycling the used product is provided by the municipal or city office.

DISPOSAL OF USED ELECTRICAL AND ELECTRONIC EQUIPMENT

Used electrical and electronic products, due to environmental protection, must not be thrown into household waste with regular municipal waste, but disposed of properly. Information on collection points and their opening hours is provided by the relevant office.

This product complies with the requirements of the relevant European and national directives that apply to it. The product meets European and national requirements on the safety of devices and products.

This product complies with the requirements of the relevant European and national directives that apply to it, RoHS.

Improper disposal of batteries/accumulators poses a risk to the environment!

Batteries/accumulators should not be disposed of with household waste. They may contain harmful heavy metals and should be treated as special waste. The chemical symbols for heavy metals are as follows: Cd = cadmium, Hg = mercury, Pb = lead.

Therefore, used batteries / accumulators should be taken to municipal hazardous waste collection points.

The built-in battery must not be disassembled for disposal. The product must be handed over in its entirety to a collection point for used electronics.

We reserve the right to make changes to the text, design and technical data of the product without notice.

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