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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow its instructions may pose a threat to life or health.

APPLICATION AND DESCRIPTION OF THE DEVICE

The vibration massager is a pleasantly vibrating device that will exercise for you!

It will help you burn excess weight and get rid of cellulite. It will also soothe muscles, tighten sagging skin and above all allow you to lose weight quickly - without diet, sacrifice and stress.

The shaping belt is definitely the most pleasant way to lose weight and shape a nice figure. It is suitable for people of any weight and any type of figure. It allows you to quickly slim your waist and fight cellulite.

The vibration system causes muscle contraction (30-50 times per second), which in turn affects the conversion of fat into energy, which is then excreted from the cells. It can be used during various daily activities, including at work or while resting.

It is great for both women and men. It is an ideal solution for people who do not have time to exercise or go to the gym. The slimming belt is your personal mini gym that allows you to exercise muscles of almost the entire body.

Nicely shaped, modeled muscles;

Slim legs, hips, stomach, waist and arms;

Cellulite removal;

A rested, relaxed and firm body.

WARNING

Danger of improper use!

Using the device for purposes other than its intended purpose and/or using the device in any other way may involve various risks.

The device must only be used for its intended purpose.

The procedures described in this instruction manual must be followed.

6 FUNCTIONS

- Slimming
- Slimming - breaks down fat tissue
- Anti-cellulite – breaks down fat deposits – cellulitis
- Strengthens muscles - tightens the muscles, making them harder and more elastic, just like when exercising at the gym
- Massage - the perfect massage after a hard day
- Relaxing

TIP

- Using the device for purposes other than its intended purpose and/or using the device in any other way may involve various risks.
- Use the device only for its intended purposes.
- The procedures described in this operating manual must be followed. Claims of any kind resulting from damage caused by improper use, unauthorised repairs, unauthorised changes or the use of unauthorised spare parts are excluded. The risk lies solely with the user.
- If you experience any skin discomfort or irritation, discontinue use immediately.

TECHNICAL DATA

- The belt works on the principle of vibration,
- Power supply - 230V,
- Remote controlled to select operating mode and intensity level;
- It has 2 operating modes: (automatic or manual)
- It has 5 intensity levels;
- Length: max. 125 cm;
- The belt should not be used in the following cases:
- By people with heart disease, coronary artery disease or with a pacemaker,
- By pregnant women.
- Additionally, do not place the device on sensitive body parts, etc.

THE SET INCLUDES:

Brand new massage belt

Controller-pilot

Charger

Measure

Bag

User manual

Original packaging

SAFETY INSTRUCTIONS

- This product may be used by children aged 8 years and over and by persons with reduced physical, sensory or mental capabilities or lack of experience or knowledge, provided they are supervised or have been instructed on the safe use of the product and understand the resulting risks.
- Never immerse the device in water.
- For cleaning, use a damp cloth or a mild detergent.
- Do not use a damaged device.
- The product is not intended for use by children.
- Children must not play with the appliance.
- Cleaning and user maintenance shall not be performed by children without supervision.
- The product itself and its accessories are not toys, please keep them out of the reach of children to avoid danger.
- Make sure that packaging materials are not left unattended. Children may start playing with them, which is dangerous.
- Protect the product from extreme temperatures, direct sunlight, strong vibrations, high humidity, moisture, flammable gases, vapors and solvents.
- Do not expose the product to mechanical stress.
- If safe operation is no longer possible, discontinue use and secure the product against further use. Safe operation is no longer possible if the product: - has been damaged, - does not function properly, - has been stored for an extended period under unfavourable conditions or - has been excessively stressed during transport.
- It is forbidden to use the product if any of its parts are damaged. In the event of damage to the cable, it is forbidden to make repairs yourself.
- Do not use if the power cord or power adapter is damaged or if it is not properly plugged into the outlet.
- Do not disassemble the device yourself.
- Do not hit or drop this product. Stop using it if it is damaged in a fall or other mechanical damage.
- Store the device in a cool and dry place.
- Before using the device, make sure that the mains voltage matches the required operating voltage of the device.
- Always use the product as intended.

ELECTRIC SHOCK HAZARD

- Connect the device only to a properly installed and grounded power outlet. The mains voltage must match the data on the nameplate of the device / mains adapter.
- Make sure that the power cable does not get wet or damp during operation. The cable must be routed so that it is not pinched or damaged.
- Keep the power cable away from hot surfaces.
- Do not perform any repairs to the device yourself.
- Any repairs may only be performed by service or by qualified specialists.
- Before cleaning the device, unplug the power supply from the wall outlet.
- Do not use the device in a damp or wet environment, such as a bathroom.
- Have a damaged plug or power cable repaired immediately by a qualified specialist or service center in order to avoid any danger.
- Do not insert any needles or sharp objects into the device.

WARNING! INJURY RISK!

- This device must not be used to massage infants, small children, helpless persons or animals.
- Never use this device if you experience dislocation or pain in the area being massaged.
- Stop the massage immediately if it becomes uncomfortable or causes pain.
- Do not use the device on swollen, burned, irritated, diseased or injured areas of the skin or body.
- If you have any doubts regarding use, please consult your doctor.
- The device must not be used on the face (e.g. eyes), larynx, feet or other sensitive parts of the body.
- Do not cover the device (pillow, blanket, etc.) during use. Do not use the device near gasoline or other flammable materials.

WARNING – MATERIAL DAMAGE!

- Use the appliance indoors only.
- Only use the device with the supplied mains adapter.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a waste disposal point designated by local authorities. Information on the possibilities of recycling the used product is provided by the municipal or city office.

DISPOSAL OF USED ELECTRICAL AND ELECTRONIC EQUIPMENT

Used electrical and electronic products, due to environmental protection, must not be thrown into household waste with regular municipal waste, but disposed of properly. Information on collection points and their opening hours is provided by the relevant office.

This product complies with the requirements of the relevant European and national directives that apply to it. The product meets European and national requirements on the safety of devices and products.

This product complies with the requirements of the relevant European and national directives that apply to it, RoHS.

Warning of electric shock! Danger to life!

We reserve the right to make changes to the text, design and technical data of the product without notice.

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