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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow them may pose a risk to life or health.

APPLICATION

The Air Compression Calf Massager is a device designed for a relaxing leg massage, based on gentle air compression technology. It works by rhythmically pressing and relaxing the muscles, which promotes blood circulation, reduces tension, and eliminates the feeling of tired and heavy legs.

This device is intended for home use only and is not intended for medical or commercial use, and is not intended to replace medical therapy. The device is intended for indoor use in dry rooms. It is not intended for medical, therapeutic, or commercial use.

For safety and CE certification reasons, this product must not be rebuilt or modified in any way. Using the product for purposes other than those described may result in damage to the product. Improper use may also result in hazards such as short circuits, fire, electric shock, etc.

WARNING

Danger of improper use!

Using the device in a manner other than its intended purpose and/or using the device in any other way may involve various risks.

The device should only be used for its intended purpose.

Please follow the procedures described in this instruction manual.

Note : It is recommended to consult a doctor before using the massager, especially if you have any problems in the calf area!

ACTIVATION

- Note: Packaging material must be completely removed from the product.
- Make sure the device components are properly assembled.
- Before first use, charge the device using the USB cable. Connect the USB cable to a USB power source and connect the USB cable plug to the product's USB input port.
- Once charged, disconnect the cable and turn on the device using the power button.
- Do not turn on the device while it is connected to the charging cable.
- Place the massager on your calf, adjusting it to your leg using the adjustable Velcro straps.
- Hold the ON/OFF button for 3 seconds to turn on the device. By default, it will turn on *SOOTHING mode*.
- To activate the heating function, press the HEATING button . Each subsequent press increases or decreases the heat level from 1–3.
- The air massage mode can be changed using the *MODE button* – three intensity levels are available: *soothing* , *tira* and *powerful* .
- Vibration function is activated by the *VIBRATION button* , also in three intensity levels.
- The device automatically turns off after 20 minutes of operation to ensure safety and energy savings.
- After finishing the massage , turn off the device and then store it in a dry place.

TIP

- Using the device in a manner other than its intended purpose and/or using the device in any other way may involve various risks.
- Use the device only for its intended purposes.
- The procedures described in this operating manual must be followed. Claims of any kind resulting from damage caused by improper use, improper repairs, unauthorized modifications, or the use of unauthorized spare parts are excluded. The risk lies solely with the user.
- If you experience any discomfort or skin irritation, discontinue use immediately.

TECHNICAL DATA

- Intensity levels: 3
- Heating control: 3 levels
- Massage mode: soothing / tira / powerful
- Power supply: battery
- Battery capacity: 2000mAh
- Power: 7.5W
- Charging time: ~4h
- Working time: up to 20 minutes per cycle
- Charging: USB cable

SAFETY INSTRUCTIONS

- This product may be used by children aged 8 years and over and by persons with reduced physical, sensory or mental capabilities or lack of experience or knowledge, provided they are supervised or have been instructed on the safe use of the product and understand the resulting risks.
- Never immerse the device in water.
- For cleaning, use a damp cloth or a mild detergent.
- Do not use a damaged device.
- The product is not intended for use by children.
- Make sure packaging materials are not left unattended. Children may play with them, which is dangerous.
- Protect the product from extreme temperatures, direct sunlight, strong vibrations, high humidity, moisture, flammable gases, vapors and solvents.
- Do not expose the product to mechanical stress.
- If safe operation is no longer possible, discontinue use and secure the product against further use. Safe operation is no longer possible if the product: - has been damaged, - is not functioning properly, - has been stored for an extended period under unfavorable conditions, or - has been subjected to excessive stress during transport.
- Do not use the product if any part is damaged. If the cord is damaged, do not attempt to repair it yourself.
- Do not disassemble the device yourself.
- Protect the product from moisture.
- Make sure not to damage the USB cable with sharp edges or hot objects. Fully unwind the USB cable before use.
- The included USB cable is only suitable for this product.
- Turn off the product immediately and remove the cable from the product if it smells burning or smoke is visible. Have the product checked by a qualified technician before using it again.
- If you're using a USB power adapter, the power outlet must be easily accessible at all times so that you can easily unplug the USB power adapter in an emergency. Please also follow the manufacturer's operating instructions.
- Before each charging process, check the USB cable and device for any damage. Never use or charge the device if you notice any damage.
- Do not use the massager on open wounds, inflammations or in case of circulatory disorders without consulting a doctor.
- Before use, ensure that the existing mains voltage complies with the required operating voltage of the device.
- Do not cover the device while it is in operation – this may cause overheating.
- The device is intended for home use only .
- Do not use while sleeping.

SAFETY INSTRUCTIONS FOR BATTERIES/ACCUMULATORS

- Keep batteries/rechargeable batteries out of the reach of children. If swallowed, seek medical attention immediately!
- Disposable batteries must not be recharged. Batteries/rechargeable batteries should not be short-circuited and/or opened. This may result in overheating, fire, or explosion.
- Never throw batteries/rechargeable batteries into fire or water.
- Batteries/rechargeable batteries should never be subjected to mechanical stress.
- Risk of acid leakage from batteries/accumulators.
- Avoid extreme conditions and temperatures that may affect batteries/accumulators, e.g. radiators/direct sunlight.
- If batteries/rechargeable batteries leak, avoid contact with skin, eyes, and mucous membranes! Immediately flush the affected areas with clean water and seek medical attention!
- Leaking or damaged batteries/rechargeable batteries can cause chemical burns when in contact with skin. Therefore, wear appropriate protective gloves.
- The product has a built-in battery that cannot be replaced by the user.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a designated waste disposal site designated by local authorities. For information on recycling options for used products, please contact your local municipality or city office.

DISPOSAL OF USED ELECTRICAL AND ELECTRONIC EQUIPMENT

For environmental reasons, used electrical and electronic products should not be disposed of with regular household waste, but should be disposed of properly. Information on collection points and their opening hours can be obtained from the relevant authorities.

This product complies with the requirements of the applicable European and national directives. The product meets European and national requirements for the safety of devices and products.

This product complies with the requirements of the relevant European and national RoHS directives.

Improper disposal of batteries/accumulators poses a threat to the environment!

Batteries/accumulators should not be disposed of with household waste. They may contain harmful heavy metals and should be treated as special waste. The chemical symbols for heavy metals are: Cd = cadmium, Hg = mercury, Pb = lead.

Therefore, used batteries/accumulators should be taken to municipal hazardous waste collection points.

The built-in battery must not be disassembled for disposal. The product should be taken to a collection point for used electronics.

We reserve the right to make changes to the text, design and technical data of the product without notice.

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