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Dear Sir or Madam, thank you for purchasing our product!

Before using the product, please read the following instructions to ensure proper use of the product.

Please keep this manual for future reference and follow its recommendations, as failure to follow its instructions may pose a threat to life or health.

APPLICATION AND DESCRIPTION OF THE DEVICE

A practical shoulder and arm band with a heating function is an ideal solution for people seeking relief from pain, muscle tension or joint ailments. Equipped with a 3-stage temperature regulation system, it allows you to adjust the heat level to individual needs, supporting regeneration and improving circulation in the arm and shoulder area.

The band is powered by a USB-C port, which makes it convenient and mobile - you can use it at home, in the office or on the go, connecting it to a power bank, laptop or adapter. The ergonomic shape and adjustable straps ensure a comfortable fit to the body, and the soft, breathable material guarantees comfort even for a long time.

This multifunctional device is a great support in everyday relaxation and in the therapy of pain resulting from injuries, overloads or tension. Ideal for athletes, physically active people, as well as anyone who wants to take care of for health and comfort. The headband combines modern technology with practical design, being a great gift idea.

COMMISSIONING/INSTALLATION INSTRUCTIONS

- Make sure the band is clean and dry before first use.
- Check that all items are complete:
- Armband with adjustable straps,
- USB power cable,
- User manual (if included).
- Find a comfortable place to put on the headband – preferably while sitting or standing in front of a mirror.
- Place the band on your arm so that the main heating part is on the area requiring heat therapy.
- Adjust the straps to ensure a proper fit – they should not be too loose or too tight
- Thread the strap through the adjustable clasps and fasten with the Velcro.
- Make sure the band fits snugly to your body but does not restrict blood circulation.
- Connecting the power source
- Connect the USB cable to the wristband.
- Connect the other end of the cable to a power bank, computer, USB charger or other 5V power source.
- Press and hold the power button for 3 seconds to turn on the device.
- The LED indicator will light up to show that the wristband is working.
- Temperature control
- Change heating modes by pressing the button briefly:
 - Red (high temperature)
 - Yellow (Medium Temperature)
 - Blue (Low Temperature)
- Press and hold the button for 3 seconds to turn off the device.
- Unplug the power cable.
- Carefully remove the band and store it in a dry place.
- Note: Do not use the band for more than 30-60 minutes in one session to avoid overheating the skin.

TIP

- The device must be checked for completeness of delivery and any visible damage.
- In the event of an incomplete delivery or damage due to faulty packaging or transport, please contact customer service.

TECHNICAL DATA

- 3 heating modes
- Power supply: USB-C (cable included)
- Cable length: 100cm
- Color: black
- MEASUREMENTS (measured flat, without stretching):
- Length: 35cm
- Additional strap length: 82cm
- Width: 35cm

SAFETY INSTRUCTIONS

- This product may be used by children aged 8 years and over and by persons with reduced physical, sensory or mental capabilities or lack of experience or knowledge, provided they are supervised or have been instructed on the safe use of the product and understand the resulting risks.
- Children are not allowed to play with the product.
- Cleaning and maintenance may not be performed by children without supervision.
- Never immerse the device in water.
- For cleaning, use a damp cloth or a mild detergent.
- Do not use a damaged device.
- Make sure that packaging materials are not left unattended. Children may start playing with them, which is dangerous.
- Protect the product from extreme temperatures, direct sunlight, strong vibrations, high humidity, moisture, flammable gases, vapors and solvents.
- Do not expose the product to mechanical stress.
- If safe operation is no longer possible, discontinue use and secure the product against further use. Safe operation is no longer possible if the product: - has been damaged, - does not function properly, - has been stored for an extended period under unfavourable conditions or - has been excessively stressed during transport.
- It is forbidden to use the product if any of its parts are damaged. In the event of damage to the cable, it is forbidden to make repairs yourself.
- Do not disassemble the device yourself.
- Protect the product from moisture.
- The light elements cannot be replaced. If the light elements have stopped working due to wear, the entire product must be replaced.
- Make sure not to damage the USB cable with sharp edges or hot objects. Unwind the USB cable completely before use.
- The included USB cable is only suitable for this product.
- Immediately turn off the product and remove the cable from the product if it smells burning or smoke is visible. Have the product checked by a qualified technician before using it again.
- If a USB power supply is used, the socket must always be easily accessible so that the USB power supply can be easily removed from the socket in an emergency. Please also observe the manufacturer's operating instructions.
- Do not look into the beam of light.
- Due to the extremely rapid speed of light change, visual disturbances may occur in children or particularly susceptible individuals (e.g. those with mental illnesses).
- Before use, make sure that the existing mains voltage complies with the required operating voltage of the device.
- Disconnect the USB cable before assembly, disassembly and cleaning.
- Always use the product as intended.

BATTERIES/ACCUMULATORS SAFETY INSTRUCTIONS

- Only connect to compatible power sources – e.g. power banks, computers or 5V USB adapters.
- Do not use damaged cables – replace the cable if you notice any abrasions or damage.
- Do not immerse the band in water or use it in a humid environment to avoid the risk of electric shock.
- Do not open or modify the device - this may void the warranty and pose a health risk.
- Do not use the bandage on skin areas with wounds, burns, irritations or inflammation.
- Do not use the band for more than 30-60 minutes in one session to avoid overheating the skin.
- Do not use while sleeping - there is a risk of overheating and excessive pressure on the arm.
- Discontinue use if you experience discomfort, pain or excessive skin heat.
- After use, unplug from power source and allow to cool before storage.
- Only clean your band with a damp cloth – do not machine wash it or submerge it in water.
- Store in a cool, dry place, away from direct sunlight and heat sources.
- Do not fold or crush the cables - this may damage them.
- Do not use the bracelet if you have a pacemaker or other medical devices that operate on electricity.
- People suffering from diabetes, circulatory disorders, skin diseases or are hypersensitive to heat should consult a physician before use.
- Do not use the band on children, the elderly or people with limited ability to feel heat without consulting a doctor.
- Do not use during pregnancy unless your doctor gives you permission.

TIPS AND INFORMATION ON THE MANAGEMENT OF USED PACKAGING

The packaging is made of environmentally friendly materials that can be disposed of at your local recycling center.

Used packaging material should be delivered to a waste disposal point designated by local authorities. Information on the possibilities of recycling the used product is provided by the municipal or city office.

DISPOSAL OF USED ELECTRICAL AND ELECTRONIC EQUIPMENT

Used electrical and electronic products, due to environmental protection, must not be thrown into household waste with regular municipal waste, but disposed of properly. Information on collection points and their opening hours is provided by the relevant office.

This product complies with the requirements of the relevant European and national directives that apply to it. The product meets European and national requirements on the safety of devices and products.

This product complies with the requirements of the relevant European and national directives that apply to it, RoHS.

Do not look into a light source when it is on. It may be harmful to your eyes.

Improper disposal of batteries/accumulators poses a risk to the environment!

Batteries/accumulators should not be disposed of with household waste. They may contain harmful heavy metals and should be treated as special waste. The chemical symbols for heavy metals are as follows: Cd = cadmium, Hg = mercury, Pb = lead.

Therefore, used batteries / accumulators should be taken to municipal hazardous waste collection points.

The built-in battery must not be disassembled for disposal. The product must be handed over in its entirety to a collection point for used electronics.

We reserve the right to make changes to the text, design and technical data of the product without notice.

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